

Det-Sol

500

HOSPITAL GRADE
BLEACH
DISINFECTANT

Suitable for general sanitation

RINSE FREE

1 SACHET MAKES 1 LITRE

Contains Sodium Dichloroisocyanurate 10-30%

6g NET

Eucalip Bio-Chemical Group Pty Ltd
www.eucalipgroup.com
Phone: 1300 880 739
Fax Direct: 1300 880 749
eucalipgroup@bigpond.com

CAUTION: Highly reactive; may cause fire or explosion

AUST L 90768



 +  = 500ppm

Directions

- Add one entire sachet to one litre of water, stir and fully dissolve to produce 500ppm available chlorine for sanitation and cleaning of surfaces, benches, walls, floors, etc
- Do not allow sachet to get damp
- Use solution within 24 hours of preparation
- Not for use on medical devices

Safety

- Do not mix with any other chemicals or detergents
- Wear suitable gloves and avoid breathing dust
- Avoid contact with skin and eyes

First Aid

- If in contact with skin or eyes, flush with large quantities of water
- If swallowed, seek immediate medical attention
- Refer to Material Safety Data Sheet on www.eucalipgroup.com

Fax Direct: 1300 880 749

5 year shelf life from DOM